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Instruction for Use

1. Indications for Use:

Defibrillation Electrodes are also called Multifunction Electrode Pads, which are accessories for defibrillator devices. It is applied on the patient's chest and will be used for administering cardiac defibrillation shock (defibrillation and cardioversion), monitoring and providing transcutaneous cardiac pacing when medically indicated. It is available for from pre-school children to Adults.

2. Warnings:

- Not following these instructions or any other misuse or misapplication of any electrode may result in patient burns or ineffective therapy. (Note: Reddening of skin is normal.)
- Use these electrodes only with compatible defibrillator devices as listed on the IFU. Use of this product with any device not listed may result in potentially hazardous connection.
- Strangulation hazard – Route cables away from patient's throat.
- Non-invasive pacing may cause skin irritation or burns. Extensive pacing may lead to air bubbles. To ensure optimal contact, smooth out bubbles by hand after applying the electrode.
- The device shall be kept well clear of other electrodes or metal parts in contact with patient.
- The operator shall be advised that other device which has no defibrillation-proof applied parts shall be disconnected from the patient during defibrillation.
- The pregnant or nursing women shall be used the device under the guidance of the medical care personnel.

3. Safety precautions:

- Do not use if product is expired.
- Do not crush, bend or fold electrodes or store them under heavy objects.
- Ensure patient's chest is clean and dry.
- Do not open pouch until ready for use.
- Inspect pouch, electrode, cable, connector and the environmental storage limitation of device prior to use.
- Do not use, if compromised.
- Use electrodes within 1 day once pouch has been opened.
- Do not use, if gel is dry.
- Keep applied defibrillation electrodes clear of any other electrodes or metal parts in contact with patient. Avoid electrode placement above an internal pacemaker.
- Do not use additional gel on electrodes.
- Do not overlap defibrillation electrodes.
- Avoid electrode placement above an internal pacemaker.

9

- Do not reposition electrodes and do not touch patient during defibrillation.
- Do not discharge hand-held paddles through these electrodes.
- Do not re-use the defibrillation electrodes. If re-used, the electrical properties might be insufficient, which could lead to a patient injury/hazard.
- Adult defibrillation electrodes are suitable for patients older than 8 years old or weighing more than 25 KG; In an emergency, do not delay treatment to determine the patient's actual weight.
- Reuse results in a risk of cross-infection from one patient to another.
- Use with CF type defibrillator monitor.
- Use separate ECG electrodes and cable if need monitoring.
- Always keep a second pack of defibrillation electrodes with the defibrillator.
- Do not use multifunction electrodes for magnetic resonance imaging (MRI) procedures.
- The operators shall avoid contact between parts of the patient's body such as exposed skin of head or limbs, conductive fluids such as gel, blood or saline and metal objects such as a bed frame or a stretcher which may provide unwanted pathways for the defibrillating current.
- For use only for patient from children under 8 years old or weighing less than 55 pounds (25 kg).
- These electrodes are single-patient use. Reuse per unit could decrease the conductive function and increase the risk of the burning.
- Replace per electrodes set after same patient attachment:
 - 50 defibrillation shocks
 - pacing reaches 8 hours in: 100mA/80bpM/40ms pulse duration;
 - 24 hours on skin
 - In critically ill patients or those under contact, airborne, or droplet precautions.

4. Instruction for Use patients:

- Select placement sites on body of patient according to defibrillator manufacturer's operation instructions for Anterior-Anterior electrode placement (1a).

- Anterior-Posterior electrode placement for pre-school children if suitable.
- Ensure skin site is clean, dry and hairless. Shave if required.
- Remove electrodes from package (3a) and then remove protective cover from electrodes (3b).
- Defibrillation electrodes are applied to the appropriate place on the patient's skin. There are two methods to attach the defibrillation electrode plate:
 - Anteroposterior (B):** One electrode plate is attached to the subscapular area of the back, and the other one is attached to the left edge of the sternum 3–4 intercostal level.

10

- Anterior lateral position (A):** One electrode plate is attached to the right edge of the sternum between the 2–3 intercostals (the base of the heart) and the other to the 5th intercostal area (apex) in the anterior line of the left axilla.

The distance between the two electrode plates should not be less than 10 cm. The anteroposterior position of defibrillation electrodes is recommended for children (B). The actual application method shall be judged by the user according to the patient's body size.

- Press down firmly to ensure good contact of the entire adhesive surface with the skin. Connect the electrodes to the defibrillator.
- Proceed with defibrillation/pacing according to defibrillator manufacturer's operating manual.
- Disconnect electrode connector from defibrillator.
- Remove electrodes by slowly peeling electrode from the patient's skin. Dispose used electrodes following local state or federal guidelines.

5. Indications for Use:

Series No.	Recommended defibrillator
A-I-HP-120-01 A-O-HP-120-01 C-I-HP-120-01 C-O-HP-120-01	HEARTSTART FR2, FR2+, FR3, MRx, XL, XL+, HEARTSTREAM FORERUNNER
A-I-QC-120-01 A-O-QC-120-01 C-I-QC-120-01 C-O-QC-120-01	PHYSIO-CONTROL LIFEPAK 9, 10, 12, 15, 20, 500, 1000; BENEHEART D3, BENE-HEART D6; Jousing
A-I-ZO-120-01 A-O-ZO-120-01 C-I-ZO-120-01 C-O-ZO-120-01	AED7000, Zoll serie M-, E-, X-, PD 1200, 1400, 1600, 1700, 2000
A-O-ZO-120-02 A-I-ZO-120-02 C-O-ZO-120-02 C-I-ZO-120-02	Zoll serie M-, E-, R-Series AED Plus, AED Pro
A-I-NK-120-01 A-O-NK-120-01 C-I-NK-120-01 C-O-NK-120-01	Nihon Kohden
A-I-SC-120-01 A-O-SC-120-01 C-I-SC-120-01 C-O-SC-120-01	Schiller: DG4000, DG5000, DG6002, Skity, OBLS300, OBLS300A
A-I-CU-120-01 A-O-CU-120-01 C-I-CU-120-01 C-O-CU-120-01	CU Medical Systems, Cmos Drake, Progetti, Cardiac Science, GE

11

6. Shelf life

Shelflife & Production date see package label.

7. Contraindications

None.

8. Side-effects

Skin redness of burning if improper usage.

9. Environment condition

- A Storage/transportation ambient temperature range of -10°C ~ +40°C
- An operation ambient temperature range of 15°C ~ 40°C
- An environment relative humidity range of ≤80%
- An environment atmospheric pressure range of 50kPa to 106kPa

12. Explanation for used symbols


Table 1

Table with 3 columns: Emissions test, Compliance, and Electromagnetic environment - guidance. It details RF emissions (CISPR 11), harmonic emissions (EN61000-3-2), and voltage fluctuations (EN61000-3-3) with their respective compliance levels and environmental requirements.

13

Table 2

Table with 4 columns: Immunity test, EN60601 test level, Compliance level, and Electromagnetic environment - guidance. It covers electrostatic discharge (ESD), electrical fast transient/burst, surge, voltage dips, and power frequency magnetic fields.

14

Table with 4 columns: Immunity test, EN60601 test level, Compliance level, and Electromagnetic environment - guidance. It details conducted and radiated RF immunity tests, including field strength and separation distance requirements.

15

Table 4

Table with 4 columns: Rated maximum output power of transmitter (W), Separation distance according to frequency of transmitter (m), and guidance. It provides separation distance requirements for different transmitter power levels and frequencies.

13. Reporting incidents involving the use of the device

Any serious incident related to the device must be reported to the manufacturer and to the authority of the Member State appropriate for the user's or patient's place of residence.

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
Table with 2 columns: Seriennummer and Empfohlener Defibrillator. Lists various device models and their compatible defibrillators.

19

Table with 2 columns: Seriennummer and Empfohlener Defibrillator. Lists additional device models and their compatible defibrillators.

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