## **PL Elektrody EKG**

**EN ECG electrodes**

**DE EKG-Elektroden**

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**SK EKG elektródy**

**CS Elektrody EKG**

**RO Electrozi ECG**

**ET EKG elektroodid**

**LV EKG elektrodi**

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EN

Instruction for disposable ECG electrodes, EK-S type.

## 1. Intended use and operation of the device

ECG electrodes are intended for diagnostic tests and monitoring of the patient's heart functions.

Electrodes are placed on the skin to transmit an electrical signal from the surface of the body via ECG cables and leadwires to the ECG apparatus or cardiac monitor, which convert the signal into electrocardiogram or vector cardiogram.

## 2. Patient target groups

Electrodes are dedicated to use by qualified medical personnel in infants, children and adults, as suggested by the manufacturer in point 12.

## 3. Contraindications

There are no known contraindications for use of the device.

## 4. Possible side effects

ECG electrode is a product that comes into direct contact with the patient's skin. There is a risk of irritation and allergic reactions in people with extremely sensitive skin.

The occurrence of skin reactions requires medical advice.

## 5. Warnings

- 5.1. Electrodes are disposable medical devices. Do not use again. Re-use may cause patients infection.
- 5.2. Use electrodes only on intact skin. Do not stick the electrode on wounded, infected or lesioned skin.
- 5.3. ECG electrodes are conditionally approved for use in a magnetic field environment, i.e., during magnetic resonance imaging (MRI), under the criteria outlined in the table in section 14. In situations that fall outside the specified criteria, leaving the electrodes on the skin may cause burns to the patient.
- 5.4. When performing electrosurgical procedures, place the electrodes as far away from the treatment site as possible to minimize the risk of burns to the patient's skin.

## 6. Cautions

- 6.1. Do not stick the electrodes on moist, lubricated skin due to the risk of the electrode coming loose.
- 6.2. Rough removal of the electrode may cause skin irritation.
- 6.3. When replacing an electrode with a new one, change the place of application of the next electrode to avoid skin irritation.
- 6.4. For disposable ECG electrodes, clamps and snaps are appropriate types of ECG cables/leadwires terminations. In case of cables/leadwires terminated with a banana plug, an adequate ECG adapter must be used.

- 6.5. Do not use electrodes that are wet or in damaged packaging. Do not use electrodes which are expired.
- 6.6. Do not open the packaging when there is no necessity of using electrodes for the tests because of the risk of the gel drying out.
- 6.7. The use of alcohol to prepare the skin for testing is not recommended due to its skin impedance increasing effect. However, if alcohol is used, make sure it has completely evaporated from the skin before applying the electrode.
- 6.8. Excessive movement of the patient adversely affects the quality of the ECG signal.
- 6.9. In case the electrode comes off of the skin, replace it with the new one.
- 6.10. Applying electrodes on excessively hairy skin may be unable for proper ECG testing.
- 6.11. Do not use abrasive agents on infants, children and people with particularly sensitive skin.
- 6.12. It is important to note that the MR Conditional labeling applies only to the electrodes and not to any other devices such as monitors, cables, or leads that may be connected to the electrodes.
- 6.13. Electrodes marked as MR Conditional may only be used with cables that are also approved for use in a magnetic field environment.
- 6.14. Electrodes and their cables must be positioned in accordance with local safety regulations for performing ECG examinations in a magnetic field environment.

## 7. Directions for use

- 7.1. Electrodes marked with „RTG” are intended for use in X-ray examination.
- 7.2. To connect the electrode with ECG apparatus, use ECG cables and leadwires compatible with equipment used.
- 7.3. Skin at the electrode site should be clean, dry and devoid from excessive hair.
- 7.4. If necessary, the skin should be shaved, washed with soap and water and dried thoroughly.
- 7.5. The oblong hole or cut in some types of electrodes is intended for attaching the leadwire.
- 7.6. Do not change the position of the electrodes once attached.
- 7.7. After the examination has ended, electrodes should be disposed according to the local regulations.
- 7.8. If a small amount of gel is left on the skin after removing the electrode, remove it with paper towel or wash it off with water.

## 8. Application

- 8.1. Open the pouch along the dotted line indicated on the packaging.
- 8.2. Determine the application sites on the patient's body. If necessary, prepare the patient's skin in accordance with point 7.3. Before sticking the electrode, always make sure that the application area is dry.
- 8.3. Connect the leadwire with electrode placed on the transparent foil.
- 8.4. After removing the electrodes from a transparent foil, apply to previously prepared skin, taking into account the information

Fig. 1. Way of leadwire and electrode connection

contained in points 5, 6 and 7. It is recommended to remove from the sticky side of the electrode the white paper element ("fingerlift"), which enables removing electrode from the foil. Apply the electrode to the skin and gently press its edges.

**8.5. For Holter electrodes:**

**8.5.1.** Without removing the electrode from the transparent foil, pull the leadwire through the hole as shown on fig. 2, starting from the sticky side towards the outer side of the electrode.

**8.5.2.** Connect the leadwire located in the hole to the electrode. Orient the snap/clamp with the leadwire upwards as shown in fig. 2. Leave approx. 3-4 centimeters of loose wire (above the hole) to prevent it from disconnecting from the electrode.

**8.5.3.** Make sure to stick the electrodes with the hole for attaching the leadwires facing downwards (fig. 2).


Fig. 2. Way of leadwire an Holter electrode connection

**9. Electrode removal**

After the examination, lift the electrode and gently, slowly detach it while holding the patient's skin with your free hand. Avoid peeling off at right angles to minimize the risk of skin irritation.

**10. Storage conditions**

**10.1.** Electrodes should be stored in their original packaging, at a temperature between 5°C and 30°C. Protect from direct sunlight, moisture and dirt.


**10.2.** The expiry date of the product is given on the packaging and applies to originally sealed products, stored in accordance with the manufacturer's recommendations.

**10.3.** Unused electrodes should be kept in the original packaging. After opening the package, the pouch should be reclosed by folding its edge.

Electrodes marked "NWSG" may be stored in the opened pouch for a maximum of 14 days.

Electrodes marked "P" or "PSG" may be stored in the opened pouch for a maximum of 30 days.

## 11. Explanation of symbols

	Suggested group of patients: adults		Indicates the date		Protect against sunlight
	Suggested group of patients: pediatric. Applies to children from infancy (2 years) to 18 years of age, depending on the patient's weight.		Production batch number		Keep dry
			Product UDI code		Temperature limit
	Liquid gel		Manufacturer		Do not re-use
	Solid gel		CE mark (conformance to the Regulation (EU) 2017/745)		Caution
	Catalogue number		Medical device		Product has an Ag/AgCl sensor
	Product index		Read the instruction for use		Product does not contain natural latex
	Pieces per packaging		Do not use if the package is damaged		Product does not contain polyvinyl chloride
	Importer		Authorized representative in the Switzerland		Product intended for use during X-ray examinations
	Product conditionally approved for use during magnetic resonance imaging (MRI).				

## 12. Suggested patient groups and intended use of ECG electrodes type EK-S

REF	Patient group	Intended use
EK-S 22 PSG	pediatric	resting examination
EK-S 25 PSG	pediatric	short-term monitoring
EK-S 30 PSG RTG	pediatric	short-term monitoring
EK-S 30 PSG	pediatric	resting examination
EK-S 35 PSG ST	pediatric, for adults	stress tests
EK-S 36 P	pediatric	long-term monitoring
EK-S 38 NWSG	pediatric	short-term monitoring
EK-S 39 PSG	pediatric	Holter
EK-S 40 NWSG	pediatric	long-term monitoring
EK-S 40 PSG	pediatric	long-term monitoring
EK-S 42 PSG	pediatric, for adults	long-term monitoring
EK-S 44 PSG	for adults	resting examination
EK-S 45 PSG RTG	for adults	short-term monitoring
EK-S 45 PSG	for adults	short-term monitoring
EK-S 48 PSG	for adults	stress tests
EK-S 50 NWSG	for adults	Holter
EK-S 50 P	for adults	long-term monitoring
EK-S 50 PSG	for adults	long-term monitoring
EK-S 52 NWSG	for adults	long-term monitoring / short-term monitoring
EK-S 55 P	for adults	short-term monitoring / stress tests
EK-S 55 PSG	for adults	long-term monitoring
EK-S 56 P	for adults	resting examination
EK-S 60 NWSG	for adults	Holter / stress tests
EK-S 60 P T	for adults	Holter
EK-S 60 P	for adults	Holter
EK-S 60 PSG XC	for adults	stress tests
EK-S 60 PSG	for adults	Holter
EK-S 61 P	for adults	long-term monitoring
EK-S 61 PSG	for adults	long-term monitoring
EK-S 62 PSG XC	for adults	stress tests

## 13. Product incident reporting

Any serious incident related to the device should be reported to the manufacturer and to the authority competent for the place of residence of the user or patient.

## 14. MRI Safety Information

### MRI Safety Information

#### MR Conditional

A patient with the SORIMEX sp. z o.o. ECG electrode EK-S may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient and/or item malfunction. Additional MR Safety information may be found at [www.sorimex.eu](http://www.sorimex.eu) or by calling +48 56 657 77 20.

Parameter	Condition of Use / Information
Static Magnetic Field Strength ( $B_0$ )	1.5 T, 3 T
Type of Nuclei	Hydrogen
Static Magnetic Field ( $B_0$ ) Orientation	Horizontal, Cylindrical Bore
Maximum Spatial Field Gradient (SFG)	40 T/m (4000 G/cm)
RF Polarization (Note: formerly called RF Excitation)	Circularly Polarized (CP)
RF Transmit Coil	Integrated Whole Body RF coil can be used. Extremities local transmit RF coils can be used since there are no ECG electrodes EK-S at the extremities local transmit RF coils level
RF Receive Coil	Any receive RF coil may be used
MR System (RF) Operating Modes or Constraints	Normal Operating Mode
Whole Body Averaged SAR	Whole Body Averaged SAR $\leq$ 2 W/kg
Scan Duration and Wait Time	Maximum 1 hour of continuous scanning with a wait 15 minutes before resuming scanning. Note: Autoscanning / Autoscan Mode is considered continuous scanning. Note: Short pauses between scan sequences are considered part of the scan time.
Item Configuration	No cables or monitoring equipment shall be connected to the electrodes during the MRI examination. The electrodes must not be used for ECG monitoring while the patient is inside the MR environment
Patient's characteristics	Several ECG electrodes EK-S may be placed on the patient, provided that a minimum distance of 3 cm is maintained between each conductive part of the electrodes (snap connector and conductive gel).
MR Image Artifact	The presence of this implant may produce an image artifact. Some manipulation of scan parameters may be needed to compensate for the artifact. The image artifact caused by the device extends approximately 34 mm from the device edge on gradient echo image, and 51 mm on spin echo images, at 1.5 T.

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