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EN ECG cables and leadwires

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1. Intended use and operation of the product

- 1.1. The ECG cables are designed to transmit electrical signals collected from a human body by ECG electrodes to a device processing those signals - an ECG unit.
Leadwires are designed to transmit electrical signals from a human body through ECG electrodes to a device processing those signals - an ECG unit, individually or through an ECG trunk cable.
- 1.2. A cable type, a way of its connection, an intended use of individual endings and leads, and other data important for correct performance of tests are specified in an operating manual of a device, for which a given cable is intended.
- 1.3. A type of a connecting wire used, its length, a type of connections on the patient's side, and external elements used are adapted to a specific type/model of a unit or a trunk cable, with which a cable is to cooperate.
- 1.4. Cables and leadwires are intended for contact with intact skin of a patient.

2. Target patient groups

ECG cables and leadwires are intended to be used by qualified medical personnel in infants, children and adults.

3. Contraindications

No known contraindications to the use of the product.

4. Possible adverse effects

No allergic reactions are observed at places where ECG cables or connecting wires are in contact with intact skin of a patient. However, such reaction cannot be completely excluded in people particularly susceptible to allergies; it should be noted that no such events have been reported to the manufacturer. Incorrect cleaning or disinfection, or a failure to perform them, as well as failure to observe hygiene procedures may cause transfer of pathogens and infections.

5. Warnings

- 5.1.** A damaged cable or leadwire must be disconnected and replaced with an intact one. ECG cables and wires must be used with units that are in a perfect operating condition. When a cable, a connecting wire or a device is damaged, immediately stop the examination. Use only in rooms with a power supply and an earthing systems in a good operating condition.
- 5.2.** Leadwires are designed for use with specific types of trunk cables or ECG units. ECG cables are designed for use with specific types of ECG units. Before using a cable or connecting wires, check them for compatibility with a device (see section 12, Compatibility).
- 5.3.** Any movement, shift and/or positioning of a patient may result in artefacts in an ECG recording, and in extreme cases, in disconnecting of a cable or wires from ECG electrodes.
- 5.4.** Strong electromagnetic fields, e.g., emitted by electrosurgical equipment, may adversely interfere with a function of an ECG cable and connecting wires. Do not use an ECG cable or wires during magnetic resonance imaging (MRI) scans.
- 5.5.** When used during defibrillation, an ECG trunk cable must be shielded against a cardiac defibrillator discharge effects. Connecting wires should be used with an ECG unit or an ECG trunk cable protected against a cardiac defibrillator discharge effects (see section 14, a pictogram "Defibrillation-proof type CF applied part"). After defibrillation, check the operational condition of a cable and connecting wires.

6. Precautions

- 6.1.** Do not use an ECG cable or connecting wires with any mechanical or electrical defects.
- 6.2.** Do not use an ECG cable or connecting wires on which any fluid was spilled.
- 6.3.** Do not modify a cable or wires in any way. Any modifications will void the warranty.
- 6.4.** Cables and wires should not be bent excessively (the bend radius should be at least 60 mm).
- 6.5.** Do not pull or yank a cable or a leadwire.
- 6.6.** A cable or a connecting wire should be routed in such way that there is no risk of its mechanical damage caused, e.g., by stepping on, driving over, placing a heavy/sharp object on, or cutting it.
- 6.7.** When connecting or disconnecting a cable, hold it by a connection (socket, plug, clamp, snap) housing. When connecting or disconnecting a wire, hold it by a connection (socket, plug, clamp, snap) housing.
- 6.8.** Protect cable and wire connections against contamination and spills.
- 6.9.** Failure to observe the above recommendations may affect safety and correct operation of the product or precision of obtained records.

7. Operating use

- 7.1.** Each time before using the product, clean and disinfect it.
- 7.2.** Connect a cable/wire plug to a socket in a compatible ECG unit or a trunk cable, switch the device on and check it for a correct operation.

- 7.3. Connect connecting wires or cable wires to ECG electrodes, and place electrodes on a patient's body in accordance with applicable medical procedures and labelling of the ECG cable.
- 7.4. After the examination, disconnect the connecting wires from the electrodes and the cable/wire plug from the device.
- 7.5. Repeat the cleaning and disinfection process.

8. Cleaning and disinfection

These activities should be performed each time before the use of a cable or wires for examination of another patient. Before cleaning or disinfection, a cable/wires should be disconnected from a device or a cable with which it operates.

8.1. Cleaning

Cable or leadwire should be cleaned by delicately wiping it with a soft cloth wetted with warm water and soap or mild detergent. Before the use, the cable/wire should be dry.

During cleaning, do not wet or immerse the cable/wire, do not wash under running water, and do not use excessive force that may break wires inside damaging the cable/wire.

Do not use any corrosive agents or other aggressive chemical formulations. Avoid any contact with strong solvents: aromatic, chlorine, ketone, ethers or esters. Do not allow any contact of cleaning agents with a terminal of a cable/wire on the electrode side and with plugs on the side of the ECG unit or the trunk cable.

8.2. Disinfection

Disinfect cables and leadwires by wiping them with formulations dedicated to disinfect plastic medical devices, on a basis of active ingredients such as isopropyl alcohol (70%), glutaraldehyde, amines, and quaternary ammonium compounds. Disinfect according to instructions of a formulation manufacturer, and leave the cable/wire to dry. Do not allow any contact of the formulation with a terminal of a cable/wire on the electrode side and with plugs on the side of the ECG unit or the trunk cable.

9. Signs of product degradation

Before each use, perform visual inspection of the product. ECG cables or connecting wires cannot be used any longer when any of the below is noticed:

- breaks in the cable jacket (external insulation) continuity, protecting conducting elements of the cable (wires, shield, connecting layer)
- deep scratches not affecting the jacket continuity, but which may cause accumulation of pathogens and prevent effective cleaning and disinfection of the product
- cracks in plastic material covering product components, making electrical connections visible
- deep scratches in plastic material covering product components that do not make electrical connections visible, but which may cause accumulation of pathogens and prevent effective cleaning and disinfection of the product
- partial or complete tearing of cables from connections/plugs/patient terminals or breaking of the jacket in those areas
- visibly bent, damaged or torn pins of a plug/socket
- damaged housing protecting the conducting elements of terminals at the patient's side, protecting against a contact with metal surfaces.

Any of the above defects disqualifies the product from further use. A product with such defects should either be repaired or disposed of.

10. Ambient conditions during operation:

- temperature: +5°C to +45°C
- relative humidity: 10% to 95%
- atmospheric pressure: 700 to 1060 hPa

11. Storage conditions

- temperature: -25°C to +55°C
- relative humidity: 10% to 85%, avoid condensation
- atmospheric pressure: 700 to 1060 hPa
- protect against sunlight, moisture and contamination

11.1. Cables and connecting wires should be stored in the specified storage conditions, to prolong the product life as far as possible.

11.2. During storage, the ECG cables or connecting wires should be loosely wound, without kinks.

12. Compatibility

A device operating together with the product is specified in a label on the product primary packaging. A table specifying ECG units compatible with the products is available at www.sorimex.eu/cablekat

13. Warranty


The product life is 10 years of the production date.

The guarantee period is 12 months from the sales date.

The general guarantee terms and conditions are available at the manufacturer website, www.sorimex.eu

The guarantee does not cover defects resulting from a failure to observe provisions of this instruction for use.

14. Description of symbols used

 REF	Product catalogue number	 index	Product index		Pieces per packaging
 UDI	Product UDI Code	 MD	Medical device	 LOT	Production batch number
	Read the instruction for use		Production date		Manufacturer
	CE mark (conformance to the Regulation (EU) 2017/745)		Caution		Marking of electrical and electronic devices according to EN 50419
	Keep dry		Protect against sunlight		Acceptable temperature limits
	Defibrillation-proof type CF applied part				

15. Notification of incidents related to the product use

Each serious incident related to the product should be notified to the manufacturer and to an authority of the Member State competent for the user or patient's place of residence.

16. Product disposal

Used cable/connecting wire should be disposed of in accordance with internal hygiene procedures in force at a given health care facility and with local regulations.

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- **EKG-Diagnose-Kabel** vom Typ: KA, W.KAB.KA, KAD, W.KAB.AD, KAE, W.KAB.AE, KAP, W.KAB.AP, KAQ, KAX, W.KAB.AX, KB, W.KAB.KB, KG, W.KAB.KG, KP, KS, W.KAB.KS, KT, W.KAB.KT
- **EKG-Monitoring-Kabel** vom Typ: KA, W.KAB.KA, KAD, W.KAB.AD, KAE, W.KAB.AE, KAK, W.KAB.AK, KAL, KAM, KAN, KAP, W.KAB.AP, KAS, W.KAB.AS, KAX, W.KAB.AX, KAY, KB, W.KAB.KB, KM, W.KAB.KM, Kmm, W.KAB.Mm, KS, W.KAB.KS, PP
- **EKG-Langzeit-Kabel** vom Typ: KM, W.KAB.KM, KH, W.KAB.KH, KS, W.KAB.KS, KB, W.KAB.KB
- **EKG-Leitungen** vom Typ PD, PG, W.KAB.PD

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